

# EZ O<sub>2</sub><sup>TM</sup> Adult Finger Pulse Oximeter

## FDA CLEARED/CE MARKED Go To New Heights with the EZ O<sub>2</sub><sup>TM</sup> Color Adult Finger Pulse Oximeter

Cardiovascular and pulmonary disease doesn't have to be a prison. Get out there and stay active while ensuring your oxygen levels and pulse are within your doctor's recommended range. Also a great tool for professional as well as amateur athletes to gauge their performance during training, especially at higher altitudes.

- Comfortable, lightweight and compact
- Quick and accurate measurements of: SpO<sub>2</sub> pulse rate and perfusion index
- Visual and audio alarms for increased patient safety
- Bright, color OLED display of measurements and plethysmograph
- Automatic display of data right side up
- Low power consumption
- ECG Viewer analysis software (optional)
- Auto shut-off
- Analysis software available

### Specifications

<b>Measurements:</b>	SPO <sub>2</sub> , Pulse Rate (PR), Perfusion Index (PI), Pulse Intensity
<b>Screen Rotation:</b>	Manual
<b>Visual And Audio Alarms For Low SPO<sub>2</sub> And Low And High Pulse Rate:</b>	Yes
<b>Color OLED Screen:</b>	Yes
<b>Automatic On/Off:</b>	Yes
<b>Measuring Range:</b>	Pulse Rate: 30-240bpm SPO <sub>2</sub> : 70-99% Perfusion Index: 0-20%
<b>Accuracy:</b>	SPO <sub>2</sub> : ±3% (during 70-99%) Pulse Rate: ±2 bpm or ±2% (whichever is greater) Perfusion Index: ±0.1% (during 0.2%-2%) ±1% (during 2-10%)
<b>Resolution:</b>	SPO <sub>2</sub> : 1% Pulse Rate: 1bpm Perfusion Index: 0.1%
<b>Alarm Range:</b>	SPO <sub>2</sub> low limit: 90% Pulse Rate low limit: 50 bpm, high limit: 120bpm
<b>Size:</b>	2.6 x 1.4 x 1.3 inches 66 x 36 x 33 mm
<b>Weight:</b>	2.1 oz; 60g (with batteries)
<b>Power:</b>	2, 1.5V AAA Batteries