

EZO₂TM Color Adult/Pediatric Finger Pulse Oximeter

Asthma Shouldn't Manage You, You Should Manage It with the EZO₂TM Color Adult Finger Pulse Oximeter

Take matters into your own hands for yourself or for your children and achieve the quality of life you both deserve. Breathe the free air without the anxiety associated with certain environments or asthma attacks.

- Comfortable, lightweight and compact
- Quick and accurate measurements of: SpO₂ pulse rate and perfusion index
- Visual and audio alarms for increased patient safety
- Bright, color OLED display of measurements and plethysmograph
- Automatic display of data right side up
- Low power consumption

Specifications

Measurements:	SPO ₂ , Pulse Rate (PR), Pulse Intensity
Screen Rotation:	Automatic (180 degrees)
Visual And Audio Alarms For Low SPO₂ And Low And High Pulse Rate:	Yes
Color OLED Screen:	Yes
SPO₂ Probe(s)	Pediatric (optional) Available in 3 sizes
Automatic On/Off:	Yes
Measuring Range:	Pulse Rate: 30-240bpm SPO ₂ : 70-99%
Accuracy:	SPO ₂ : ±3% (during 70-99%) Pulse Rate: ±2 bpm or ±2% (whichever is greater)
Resolution:	SPO ₂ : 1% Pulse Rate: 1bpm Perfusion Index: 0.1%
Alarm Range	SPO ₂ low limit: 90% Pulse Rate low limit: 50 bpm, high limit: 120bpm
Size:	2.6 x 1.4 x 1.3 inches 66 x 36 x 33 mm
Weight:	2.1 oz; 60g (with batteries)
Power:	2, 1.5V AAA Batteries