

# EZ O<sub>2</sub><sup>TM</sup> Wrist Pulse Oximeter

## Don't Guess. Be In The Know with the EZ O<sub>2</sub><sup>TM</sup> Color Wrist Pulse Oximeter

Knowing where you stand gives you the peace of mind to live a normal life and saves you the money and further stress associated with an unnecessary emergency room visit.

- Up to 72 hours of continuous monitoring of SpO<sub>2</sub>, pulse rate and perfusion index
- Comfortable and lightweight; ideal for patients suffering from sleep apnea
- Stored data may be transferred to a personal computer for detailed analysis and trending reports
- Visual and audio alarms for increased patient safety
- Bright, color OLED display of measurements and plethysmograph
- Low power consumption - over 50 hours of use

### Specifications

<b>Measurements:</b>	SpO <sub>2</sub> , Pulse Rate (PR), Perfusion Index (PI),
<b>Visual And Audio Alarms For Low SpO<sub>2</sub> And Low And High Pulse Rate:</b>	Yes (Programmable)
<b>Color OLED Screen:</b>	Yes
<b>Programmable Record Intervals:</b>	1 sec/1 hour, 4 sec/4 hours, 8 sec/8 hours
<b>Trending (SpO<sub>2</sub> and Pulse Rate):</b>	Yes
<b>Automatic On/Off:</b>	Yes
<b>Measuring Range:</b>	Pulse Rate: 25-250 bpm SpO <sub>2</sub> : 35-99%
<b>Accuracy:</b>	SpO <sub>2</sub> : ±2% (during 75-99%) Pulse Rate: ±2 bpm or ±2% (whichever is greater) Perfusion Index: ±0.1% (during 0.2%-2%) ±1% (during 2-10%)
<b>Resolution:</b>	SpO <sub>2</sub> : 1% Pulse Rate: 1bpm Perfusion Index: 0.1%
<b>Alarm Range</b>	All limits are adjustable. The default limits are 90% for SpO <sub>2</sub> ; 50 and 120bpm for PR
<b>Size:</b>	2.3 x 1.9 x .9 inches 59 x 49 x 22 mm
<b>Weight:</b>	2.1 oz; 60g (with batteries)
<b>Power:</b>	2, 1.5V AAA Batteries